ASHBY PRECISION ENGINEERING DRAYTON LTD

Unit 31, Nuffield Centrum Nuffield Way Abingdon Oxon OX14 1RL

Telephon Number: 01235 531279 Fax Number: 01235 535801

Component Quality Check For Part:	DSD0077-33	LLT: 7038
Item Description: 117 600		2.5m
Mary and the second		
Item Revision: 07		

TEST RESULTS	Initials	Data
Leak Test		28.3.17
Pumped and Baked To OIS LLT Syphon Pump / Bake / Test Work Instru	uction	11
Motor Check If Applicable	/	_
Flow Check @ 4.2k	55LPM	293.17
Cleanliness Check	R	-11-

	Ashby Pre-Delivery		Oi Pre-use	
Checks	Initial	Date	Initial	Date
Is The Valve + Needle Correctly Aligned On The Elbow	6	993.17		
Are The Legs Straight	1	-/-		
Are The Legs Correctly Aligned	14			
Is There Any Impact Damage To The Legs Or The Flex	1	_,_		
Is The Backing Ring Square				
is The Needle Valve Running Smoothly	A			
Are The Valve Covers in Place				
Is The Shipping Box In Good Condition	R			¥.

Checks To Be Completed Upon Fadure Of the Syphon.	Ol Post-Use		Ashby Post-Use	
	Initial	Date	Initial	Date
Is The Valve + Needle Correctly Aligned On The Elbow				
Are The Legs Straight				
Are The Legs Correctly Aligned				
Is There Any Impact Damage To The Legs Or The Flex				
ls The Backing Ring Square				
Is The Needle Valve Running Smoothly				
Are The Valve Covers in Place				
Is The Shipping Box In Good Condition				